


# easyFlow<sup>®</sup> QC Combo 1

## Hygiene- and Quality Control - made easy!

The easyFlow® QC Combo 1 kit is a simple, compact rapid test to detect viable microorganisms from different matrices cosmetic products or colors). It is fast and simple as well as robust and reliable.

You define the product volume to test: this



## The AMODIA easyFlow® QC Combo 1 kit is:

- fast:
  - speeds up microbial control and thus reduces quarantine and storage times
- unique: gives a stable and documentable result on a lateral-flow strip
- scalable: culture volume not limited

# simple:

allows a quality control at the production site without trained personnel

- reliable:
  - is based on detecting nucleic acids
- future-ready: is prepared for increasing regulatory demands

## In five simple steps from sample to result:











#### Day before

Prepare enrichment culture with volume of choice

#### 20 Minutes

For each of the three detections: Centrifuge 1 ml culture, dissolve pellet in buffer, dispense reaction reagents, add 10 µl of the suspension

#### 25 Minutes

Place reaction tubes into the thermal cycler (or incubator), start program

#### 1 Minute

Label LFDs, add hybridisation reaction

#### 15 Minutes

Dip LFDs into chromatographic buffer, read-out result

Art-No.	Product	Price
AEF-QC1	AMODIA easyFlow <sup>®</sup> QC Combo 1 Kit (25 rx.) contains: Suspending Buffer, Hybridisation Buffer, three Probe-Mixes for the Escherichia/Shigella group, <i>Pseudomonas aeruginosa</i> and <i>Staphylococcus aureus</i> , as well as Lateral-Flow Dipsticks, Chromatographic Buffer, Instruction for Use, Short Protocol	196 EUR

## If you are looking for a test that ...

is applicable for many different sample matrices

speeds up microbial control

reduces quarantie and storage times

simplifies quality control at production site

may be performed without preceeding knowledge

fits easily into existing processes

gives a permanent, easy-to-read result

needs little equipment

has no maintenance costs

detects only the specified bacteria (according to sequence analysis)


detects only viable bacteria


is based on detecting nucleic acids

has a short application time









# ... then you should test our easyFlow® QC Combo 1

Supplier: AMODIA Bioservice GmbH

Order-No.: AEF-QC1

Package size: 25 reactions

Delivery: ex-works from Braunschweig, Germany


Storage: 4 - 8°C

Stability: min. 60 days after opening

Sample material: Tested yet (a.o.): shampoo, showering gel, tooth paste, skin cream

Note: Every sample matrix has to be validated individually!

## For further information or questions, please get in contact!



**AMODIA Bioservice GmbH** Rebenring 31

38106 Braunschweig

- GERMANY -

**\*** +49 (0) 531-260 17 64

**4** +49 (0) 531-260 17 66

http://www.amodia.com